


# ProNOVA™

*Sirolimus Eluting Stent*





Manufactured By:

Vascular Concepts Limited  
2 Deans Hall Business Park  
Oak Road, Little Maplestead  
Halstead, Essex C09 2RT  
United Kingdom  
Tel: 01787 477199  
Fax: 01787 477510

Vascular Concepts Limited  
No. 19, SV Complex (1st Floor)  
Bellary Road, Hebbal  
Bangalore 560 024  
INDIA  
Tel: (080) 23438146-49  
Fax: (080) 23439205

## ProNova Technical Data

BALLOON [mm]	2.50	2.75	3.00	3.25	3.50	4.00
PRESSURE [Mpa/bar]						
0.2/2	2.42	2.66	2.92	3.17	3.42	3.92
0.3/3	2.44	2.68	2.94	3.18	3.44	3.94
0.4/4	2.46	2.71	2.96	3.21	3.46	3.96
0.5/5	2.48	2.73	2.98	3.23	3.48	3.98
0.6/6	2.50	2.75	3.00	3.25	3.50	4.00
0.7/7	2.51	2.77	3.02	3.26	3.52	4.03
0.8/8	2.52	2.79	3.06	3.27	3.56	4.07
0.9/9	2.54	2.81	3.08	3.29	3.59	4.09
1.0/10	2.56	2.84	3.11	3.31	3.62	4.11
1.1/11	2.59	2.86	3.13	3.33	3.65	4.13
1.2/12	2.61	2.88	3.16	3.35	3.68	4.17
1.3/13	2.64	2.90	3.19	3.37	3.68	4.22
1.4/14	2.66	2.92	3.22	3.40	3.74	4.26
1.5/15	2.69	2.94	3.26	3.43	3.74	4.29
1.6/16	2.72	2.96	3.29	3.45	3.80	4.34
1.7/17	2.75	2.99	3.33	3.47	3.83	4.39
1.8/18	2.79	3.00	3.36	3.49		
1.9/19	2.81	3.01	3.38			
2.0/20	2.83					

Nominal pressure at 6 BAR [0.6 Mpa]

RATED BURST PRESSURE

## ProNova Ordering Information

STENT LENGTH	STENT DIAMETER							
	2.25 mm	2.50 mm	2.75 mm	3.00 mm	3.25 mm	3.50 mm	4.00 mm	
13mm	PN22513	PN2513	PN27513	PN3013	PN32513	PN3513	PN4013	
18mm	PN22518	PN2518	PN27518	PN3018	PN32518	PN3518	PN4018	
23mm	PN22523	PN2523	PN27523	PN3023	PN32523	PN3523	PN4023	
28mm	PN22528	PN2528	PN27528	PN3028	PN32528	PN3528	PN4028	
33mm	PN22533	PN2533	PN27533	PN3033	PN32533	PN3533	PN4033	
38mm	PN22538	PN2538	PN27538	PN3038	PN32538	PN3538	PN4038	